

Mouse&Rat Tail Vein Simulator

Greatly
Contributes to
“Replacement” and
“Refinement” of
Animal Use
Alternatives



HOW TO USE
MOVIE

For Injection and Blood Collection Training Prior to Using Live Animals


Feature 1

Accurately Produced for Realistic Training Experience

- Developed with professors and professional technicians to reproduce an actual animal.
- To be used for training before performing procedures on live animals.
- Recommended for students and beginners.

Feature 2

Vein Injection Training




- Stabilizing simulator in a restraint device (pictured above).
- Insert needles into each end of tubes (blood vessels) of the tail and inject fluid. Injection is complete when the fluid comes out from the tip of the tubes.

【Development Partner】

The simulator was developed using 3D printer and soft tissue reproduction technology under academic guidance by Dr. Masaru Kawakami of Yamagata University Faculty of Engineering and Mechanical Systems.

Feature 3

Blood Collection Training



- Stabilizing simulator in a restraint device (pictured above).
- After filling tubes with simulated blood using a syringe, insert needle into either one of the tubes (blood vessels) and apply pressure to the simulated blood in the tube with the syringe. This reproduces the blood collection process.

Mice & Rat Restraint Device


Animals
enter on their
own through the

ONE WAY
SYSTEM

Patent Pending




Rat Restraint



Mice Restraint

Greatly Contributes to
“Refinement”
of Animal Use
Alternatives



Bell Shaped Entrance

Designed for smooth
entry that prevents
injuries and stress

- Minimizes animals' distress, and significantly improves efficiency
- Ideal for procedures on the tail
- Safe to sterilize in autoclaves

Product	Animal	Inner Diameter and Base Dimension (mm)
Restraint ϕ 28	Mice	Inner Diameter 28 \times L110, Base : W136 \times D100
Restraint ϕ 32	Mice	Inner Diameter 32 \times L110, Base : W136 \times D100
Restraint ϕ 35	Mice	Inner Diameter 35 \times L126, Base : W145 \times D100
Restraint ϕ 47	Rats	Inner Diameter 47 \times L176, Base : W183 \times D100
Restraint ϕ 64	Rats	Inner Diameter 64 \times L215, Base : W243 \times D126
Restraint ϕ 70	Rats	Inner Diameter 70 \times L215, Base : W243 \times D126

Winged Infusion Sets for Laboratory Animals 25G

PAT.No.6842102



Description
Video


Easy to Handle
Gentle on Animals

Feature 1

For Both Injection and Blood Collection



Blood Collection



Injection

- Connects to hematocrit tubes or syringes and are adaptable for use in various treatments.
- Connects to both plastic and glass hematocrit tubes of various widths.

Feature 2 Easily Enters Blood Vessels



Blood Collection

Injection

- Wings allow stable grip.
- Inserts approximately parallel to the animals' tail without being angled.
- Higher rate of successful venipuncture that lead to improvement of experiment precision.

Feature 3 Minimizes Damage

- As needles are made to be non-invasive on blood vessels and on the body, injection and blood collection can be performed frequently.
- The closed system protects the safety of animals by enabling direct blood collection from the caudal vein.
- Blood is collected through the natural flow of capillaries which keeps damage to red blood cells minimal and prevents hemolysis.