


# ***Non-Invasive Blood Pressure in Mice and Rats***

## **CODA™ Monitor**



- ***Systolic & Diastolic BP***
- ***Mean Blood Pressure***
- ***Heart Rate***

**\$3,450**

- Built-in computer display
- Automatically measure blood pressure in real-time
- 99% accurate to direct blood pressure & telemetry
- Clinically validated Volume Pressure Recording (VPR)

# Choose the system that's right for you!

## CODA™ MONITOR 1 ANIMAL



### Includes:

- 1 CODA Controller
- 1 Cuff Kit
- 2 Warming Pads

### Options:

- Animal Holder
- Database Software
- Laptop Computer
- On site In-Service
- Extended Warranty

## CODA™ STANDARD 1 TO 2 ANIMALS



### Includes:

- 1 CODA Controller
- 1 Controller Stand
- 1 or 2 Cuff Kits
- 1 or 2 Animal Holders
- 1 Warming Platform
- 1 Database Software

### Options:

- Additional Channel
- Laptop Computer
- On site In-Service
- Extended Warranty

## CODA™ HIGH THROUGHPUT 4, 6, OR 8 ANIMALS NETWORK UP TO 48 ANIMALS



### Includes:

- 1 CODA Controller
- 1 Controller Stand
- 4, 6, or 8 Cuff Kits
- 4, 6, or 8 Animal Holders
- 1 or 2 Warming Platforms
- 1 Database Software
- 1 Laptop Computer
- 1 On site In-Service

### Options:

- Additional Channels
- Extended Warranty

### Features

- MRI Compatible
- Measure 8 gram mice to 950 gram rats
- Diastolic BP measured, not calculated
- Easily measure dark-skinned, C57/BL6 mice
- Measure both awake and anesthetized animals

### Volume Pressure Recording

- Systolic Blood Pressure
- Diastolic Blood Pressure
- Mean Blood Pressure
- Heart Rate
- Tail Blood Volume
- Tail Blood Flow

### References

- 2008 Validation Study
- Over 400 + Peer-Reviewed Research Papers
- Published White Paper
- Video Training Protocol DVD
- Clinical Scientist Support

**Toll-Free: 888-572-8887 or 860-626-1172**

**Email: [sales@kentscientific.com](mailto:sales@kentscientific.com)**