

## Suture Material

### NYLON MONOFILAMENT

Nylon Monofilament is made up of black and clear monofilaments that are smooth and nonporous. This material has excellent tensile strength and ease of passage through tissue while reducing trauma to surrounding tissue and preventing tissue in growth.

### SILK BLACK BRAIDED


Silk Black Braided suture has optimum handling characteristics and knot security without excess capillarity.

### GUT PLAIN-ABSORBABLE

Gut Plain is manufactured from the highest quality raw material to ensure ease of handling, dependable absorption, and high tensile strength. Gut Plain absorbs within 70 days.

### GUT CHROMIC-ABSORBABLE

Gut Chromic is Gut Plain with Chromium Salt added to allow for greater absorption time. Gut Chromic absorbs within 90 days.



## Needle Types

### REVERSE CUTTING NEEDLE

Reverse Cutting Needles feature a cutting edge on the outer curve, which is beneficial when cutting through tough tissue.



### CUT TAPER NEEDLE

Cut Taper Needles feature a cutting tip and taper body. The Cut Taper Needle is good for cutting through thick tissue.



### STRAIGHT CUTTING NEEDLE


The Straight Cutting Needles are designed for smooth penetration through tissue.



### REVERSE CUTTING $\frac{3}{8}$ CIRCLE NEEDLES



### REVERSE CUTTING $\frac{1}{2}$ CIRCLE NEEDLES



### TAPER POINT $\frac{1}{2}$ CIRCLE NEEDLES

