

# GenieTouch™ Syringe Pump

Infusion/Withdrawal Dual Syringe Pump

## Standard Features

Infusion/Withdrawal  
Multi-Syringe Capability  
Selectable Functions:

- Constant Rate
  - Time Dispensing
  - Volume Dispensing
- Ramp
- Pulse
- Stepped

Automatic Positioning/Stalling Detection  
Touchscreen Control-Glove Friendly  
USB Computer Interface  
Greater than 20 kg Linear Force  
Push-Pull Capability  
Bi-directional Left or Right Operation  
2-Year 100% Parts & Labor Warranty



Order your GenieTouch™ today:

**888.572.8887**

[www.kentscientific.com](http://www.kentscientific.com)




Sci  
Chem  
Eng  
R&D  
Bio





**Kent Scientific**  
CORPORATION


- Easy-to-use touchscreen interface
- Infusion & withdrawal
- Multi-syringe support



Infusion and Withdrawal



Infusion Volume Indicator



Carriage Position Detection

Air-Tite All Plastic HWS Norm-Ject	Becton Dickinson Interm, WW, Plastipak
Becton Dickinson Glass - all types	Hamilton 1000-Series Gastight
Cadence Science MICRO-MATE glass	Ranfac
Scientific Glass Eng SGE	Sherwood Monoject Plastic
Terumo	Unimetrics Series 9000

Wide Range of Syringes

### Specifications

# of Syringes	Up to 4
Syringe Size	0.5 ul to 60 ml
User Interface	Touchscreen
Linear Force	Greater than 20 kg
Flow Rates	Min: 0.0001 ul/hr Max: 220.82 ml/min
Selectable Functions	Constant Rate Time Dispensing Volume Dispensing Ramp Pulsed Stepped
Programmable Fill/Empty	Yes
Pusher Block Stall Detection	Yes
Computer Interface	USB
Audible Alarm Indication	Yes
Carriage Position Detection	+/-0.2 mm
Warranty	2 Years



GenieTouch

Order your GenieTouch™ today:

**888.572.8887**

[www.kentscientific.com](http://www.kentscientific.com)