

# Mouse and Rat Intubation Kits



**Starting at \$675**

**Features & Benefits:**

- **Integrated anesthesia facemask**  
*Deliver isoflurane during intubation*
- **Neck positioner**  
*Align airway structures*
- **Fiber optic lighting**  
*Illuminate airway for enhanced visualization*

Order your Intubation Kit:

**888.572.8887**



**Kent Scientific**  
CORPORATION

The intubation of small animals continues to be a challenge for researchers. With an integrated anesthesia mask, the Kent Scientific intubation kit allows the delivery of isoflurane during intubation, making it easier to intubate small animals.

The ability to adjust the angle of the neck allows alignment of airway structures to visualize the vocal chords.

The fiber optic guide light illuminates the airway to enhance visualization.

## Intubation kit includes:

- Intubation platform, positions animal at 60°
- Intubation stand with adjustable head positioning
- Anesthesia mask for anesthesia delivery during intubation
- Anesthesia mask tubing & connectors (36" length)
- Fiber optic light and charger
- Intubation catheters (5-pack)
- Suture, 1-meter length
- Small spatula to position animal's tongue
- Small mirror to confirm correct placement
- Power supply

## 1-Year 100% Parts & Labor Warranty



## Integrate the SurgiSuite into your surgical set-up

*Shown with optional FIR warming pad and controller, LED lighting and magnifier kit*



The Intubation Tube Safety Wedge is designed to assist in insertion of the catheter into the trachea and minimizes the chances of inserting the cannula too deeply. During intubation, the wedge will guide the catheter tip to the proper location between the larynx and carina. Once placed, the tight fit ensures bilateral ventilation and reduces air leaks from around the catheter.



Order your Intubation Kit:

**888.572.8887**

kentscientific.com



**Kent Scientific**  
CORPORATION