



RoVent[®] and RoVent[®] Jr.

Automated ventilator



RoVent and RoVent Jr. small animal automated ventilator are ideal for mice, rats and small animals up to 1,250 grams. Both are fully automated — simply enter the animal's weight and press run.

The Rovent includes the RightTemp temperature monitor and homeothermic warming module.

Features & Benefits:

- **Fully automated with touchscreen control**
Safely ventilate using auto-calculated settings based on animal weight
- **Volume control and pressure control**
Reduce footprint with one, dual control ventilator
- **Multi-animal capability***
Greater efficiency

**Our clinical department will be pleased to review your study protocol to confirm that our product will meet your specific study requirements.*

1-Year 100% Parts & Labor Warranty (3-Year extended warranty available)

Order your ventilator:

888.572.8887



Kent Scientific
CORPORATION

Choose the ventilator that fits your research needs:

Comparison Chart RoVent® & RoVent® Jr.



	RoVent Jr. (\$1,895)	RoVent (\$3,975)
Volume ventilation	●	●
Pressure ventilation	●	●
Pulse ventilation mode	○	●
PEEP	○	●
Assist mode	○	●
Sigh breath	○	●
I:E ratio	●	●
Tidal volume – measured	●	●
Auto set-up	Limited	●
User profiles	○	●
Alarms (parameter)	○	●
Trigger in	○	●
Analog out	○	●
Digital out	○	●
Sync out	○	●
ADI integration	○	●
Data output	○	●
Internal data logging	○	●
Realtime clock	○	●
Temperature monitor and homeothermic control module	○	Included (valued at \$1,000)

Order your ventilator:

888.572.8887
kentscientific.com



Kent Scientific
CORPORATION