

PhysioSuite® with
RightTemp®

Temperature monitor and
homeothermic control module




Includes 8" x 10" pad with
two temperature probes

Only \$1,073

Features & Benefits:

- **Far infrared warming**
Penetrates deep within the body to safely warm the animal
- **Automatic control**
Self-regulate your animal's temperature
- **Safety alarms**
Respond quickly to changes in temperature



Order your RightTemp:

888.572.8887



Kent Scientific
CORPORATION

The RightTemp temperature monitor and homeothermic control module will accurately monitor and control animal temperature while safely warming your animal with far infrared heat.

PhysioSuite with RightTemp Module includes:

- Controller
- Far infrared warming pad (8"x10")
- Two temperature sensors
- Disposable sleeve protectors (10)

1-Year 100% Parts & Labor Warranty
(3-Year extended warranty available)

Product Specifications:

| | |
|------------------------|---|
| Resistance: | 2,252 ohms |
| Tolerance: | +/- 0.2° C |
| Accuracy: | 0-70° C |
| Operating Temperature: | -60° C to 150° C |
| Sensor Diameter: | Ball tip dia: 0.077" Shaft dia: 0.062" |
| Time Constant: | 15.0 seconds |
| Heat Dissipation: | 0.75 m/W (°C) |

Optional Add-On Modules:



MouseSTAT® Pulse Oximeter and Heart Rate Monitor Module*

Measure heart rates up to 900 bpm

CapnoScan® CO₂ Monitor Module

Collect accurate end-tidal CO₂ measurements

\$4,373 per module



Far Infrared Warming Pads

Infrared heat penetrates deeper than traditional warming methods, warming your animal quickly and safely.

Pads available in three sizes.

Starting at \$72



SurgiSuite Surgical Platform

The PhysioSuite with RightTemp easily integrates with the SurgiSuite surgery platform.

Warm your animal while measuring vital signs during surgery.

\$325

* The RightTemp module is included with all PhysioSuite units — a value of \$1,000

Order your RightTemp:

888.572.8887

www.kentscientific.com



Kent Scientific
CORPORATION

01.18.1