



Chemical-Resistant Base

With a stainless steel surgical surface and polycarbonate overlay for added protection

Rechargeable Far Infrared Animal Warming Pad

Recessed in center of platform for a level surgical field

Bright LED Lighting*

Connects to multi-port USB power strip

*SurgiSuite shown with optional add-ons
Exact lighting may not be as shown due to product improvements

Talk with a specialist today:
888.572.8887



Kent Scientific
CORPORATION

SurgiSuite Includes:

- (1) Easy to clean, chemical resistant base
- (1) Stainless steel surgical surface
- (1) Integrated, rechargeable far infrared warming pad
- (2) Magnetic nose cone/intubation tubing stabilizers, (1-small and 1-large)
- (4) Magnetic limb positioners (used with elastic wire)
- (1) Elastic wire positioner (1-meter)
- (1) Detachable procedure bowl
- (5) Replacement surgical field covers



Specifications:

	SurgiSuite	SurgiSuite-LG
Base dimensions:	12.90 in x 16.80 in (32.77 cm x 42.67 cm)	15.75 in x 17.80 in (40.01 cm x 45.21 cm)
Stainless steel area:	11.90 in x 7.90 in (30.23 cm x 20.07 cm)	12.75 in x 10.75 in (32.39 cm x 27.31 cm)
Warming surface area:	5.15 in x 4 in (13.08 cm x 10.16 cm)	8.75 in x 6.75 in (22.23 cm x 17.15 cm)

Optional Add-Ons:



Monitor and regulate animal temperature

The PhysioSuite® homeothermic warming and temperature control easily connects to the included warming pad.

Ask about integrating pulse oximetry.



Customize your lighting solution*

LED lighting and magnification kit plugs directly into the included USB power strip.

**Exact lighting may not be as shown due to product improvements*



Maximum control with minimum clutter

Retractor wires lock into the included fixators to allow for linear and rotational adjustment.

Order your SurgiSuite today:

888.572.8887

www.kentscientific.com



Kent Scientific
CORPORATION