


PhysioSuite® with

# RightTemp®

Temperature monitor and  
homeothermic control module



Includes 8" x 10" pad with  
two temperature probes



### Features & Benefits:

- **Far infrared warming**  
*Penetrates deep within the body to safely warm the animal*
- **Automatic control**  
*Self-regulate your animal's temperature*
- **Safety alarms**  
*Respond quickly to changes in temperature*

Order your RightTemp:

# 888.572.8887



**Kent Scientific**  
CORPORATION

The RightTemp temperature monitor and homeothermic control module will accurately monitor and control animal temperature while safely warming your animal with far infrared heat.

## PhysioSuite with RightTemp Module includes:-


- Controller
- Far infrared warming pad (8"x10")
- Two temperature sensors
- Disposable sleeve protectors (10)

**1-Year Parts & Labor Warranty**

## Product Specifications:

Resistance:	2,252 ohms
Tolerance:	+/- 0.2° C
Accuracy:	0-70° C
Operating Temperature:	-60° C to 150° C
Sensor Diameter:	Ball tip dia: 0.077" Shaft dia: 0.062"
Time Constant:	15.0 seconds
Heat Dissipation:	0.75 m/W (°C)

## Optional Add-On Modules:



### MouseSTAT® Pulse Oximeter and Heart Rate Monitor Module\*

Measure heart rates up to 900 bpm

### CapnoScan® CO2 Monitor Module

Collect accurate end-tidal CO2 measurements



### Far Infrared Warming Pads

Infrared heat penetrates deeper than traditional warming methods, warming your animal quickly and safely.

Pads available in two sizes.



### SurgiSuite Surgical Platform

The SurgiSuite surgery platform easily integrates with PhysioSuite or any of our far infrared warming controllers.

Warm your animal while measuring vital signs during surgery.

**\* The RightTemp module is included with all PhysioSuite units - a value of \$1,000**

Order your RightTemp:  
**888.572.8887**  
kentscientific.com



**Kent Scientific**  
CORPORATION

09.22.1.NP