



RoVent[®] and RoVent[®] Jr.

Automated ventilator



RoVent and RoVent Jr. small animal automated ventilator are ideal for mice, rats and small animals up to 1,250 grams. Both are fully automated — simply enter the animal's weight and press run.

The Rovent includes the RightTemp temperature monitor and homeothermic warming module.

Features & Benefits:

- Fully automated with touchscreen control
Safely ventilate using auto-calculated settings based on animal weight
- Volume control and pressure control
Reduce footprint with one, dual control ventilator
- Multi-animal capability*
Greater efficiency

**Our clinical department will be pleased to review your study protocol to confirm that our product will meet your specific study requirements.*

Standard 1-Year Parts & Labor Warranty (Premium 3-year warranty available).
Visit our website for more information.

Order your ventilator:
888.572.8887



Kent Scientific
CORPORATION

Choose the ventilator that fits your research needs:

Comparison Chart RoVent® & RoVent® Jr.



	RoVent Jr.	RoVent
Volume ventilation	●	●
Pressure ventilation	●	●
Pulse ventilation mode	○	●
PEEP	○	●
Assist mode	○	●
Sigh breath	○	●
I:E ratio	●	●
Tidal volume — measured	●	●
Auto set-up	Limited	●
User profiles	○	●
Alarms (parameter)	○	●
Trigger in	○	●
Analog out	○	●
Digital out	○	●
Sync out	○	●
ADI integration	○	●
Data output	○	●
Internal data logging	○	●
Realtime clock	○	●
Temperature monitor and homeothermic control module	○	Included

Order your ventilator:
888.572.8887
kentscientific.com



05.23.1.NP